

<sup>1</sup>University of Miami Miller School of Medicine; Cognivue Advisory Board Member; <sup>2</sup>University of Miami Miller School of Medicine; <sup>3</sup>Cognivue Working Group, Cognivue Inc., Victor, NY. \*Presenting author.

## KEY TAKEAWAY: Cognivue Clarity<sup>®</sup> discriminated healthy controls from those with AD biomarkers, and had strong psychometric properties.

# BACKGROUND

Alzheimer's Disease (AD) and Mild Cognitive Impairment (MCI) treatment and prevention trials are underway; however, recognition of potentially eligible individuals can be challenging, which leads to delayed recruitment and high screen failure. Digital biomarkers that can accurately identify MCI and AD and correspond to AD biomarkers could greatly facilitate enrollment and track disease progression and/or response to treatment.

Cognivue *Clarity*<sup>®</sup> is an FDA-cleared 10-minute, digital cognitive testing platform using adaptive psychophysics to capture global and domain-specific performance (Figure 1). The test generates a report that can be used by clinicians to assess cognitive status.

The Bio-Hermes Study generated a firstof-its-kind, "head-to-head" comparison of a broad set of AD biomarkers. The Cognivue Clarity device was included as part of the study's assessments to explore the utility of several leading biomarkers.



Figure 1. The Cognivue Clarity device.

# **DESIGN/METHODS**

The Bio-Hermes study, funded by the Global Alzheimer's Platform Foundation<sup>®</sup>, enrolled healthy controls, as well as patients with MCI and AD. Assessments included sociodemographics, Cognivue Clarity, plasma AD biomarkers, and amyloid positron emission tomography (PET) scans. Cognivue Clarity's psychometric properties, diagnosis discrimination, correlations with AD biomarkers, and diagnostic power were analyzed.

# RESULTS

The entire Bio-Hermes study included 1,001 participants, including 417 controls, 311 participants with MCI, and 273 participants with probable AD. Participants had a mean age of 72.0.9+6.7 years and 15.4±2.8 years of education. 56.2% of participants were female and 23.8% were from underrepresented minorities (75.6% Non-Hispanic White, 11.3% Black, 11.4% Hispanic, 1.8% Asian).

The mean Mini–Mental State Examination (MMSE) score was 26.6+2.9 and mean Cognivue *Clarity* score was 62.9+17.4. Cronbach's alpha among the whole sample was 0.879 (95% CI: 0.867-0.890), supporting high reliability.

Cognivue Clarity scores differed between ethnoracial groups in the study (Table 1).

# **Cognivue Clarity<sup>®</sup> in the Detection of Biomarker-Confirmed Mild Cognitive Impairment and Alzheimer's Disease**

James E. Galvin, MD, MPH<sup>1</sup>; Michael J. Kleiman, PhD<sup>2</sup>; Heather M. Harris<sup>3,\*</sup>; Paul W. Estes<sup>3</sup>

	Non-Hispanic White	Non-Hispanic Black	Hispanic	p-value		Non-Hispanic White	Non-Hispanic Black	Hispanic	p-valu
Age	72.5 (6.6)	70.2 (6.2)	71.2 (6.9)	.001ª	Age	72.3 (6.7)	69.9 (6.2)	70.7 (6.2)	.042ª
Cognivue	63.9 (17.4)	59.1 (14.4)	58.2 (18.6)	<.001 <sup>b</sup>	Cognivue	64.8 (17.5)	62.3 (13.6)	59.6 (18.4)	.195
MMSE	27.0 (2.7)	25.3 (3.0)	25.2 (3.2)	<.001 <sup>b</sup>	MMSE	26.9 (2.9)	25.8 (3.1)	25.9 (3.5)	<.001°
Centiloid Level	30.2 (47.9)	23.4 (37.9)	29.9 (45.8)	.399	Post-hoc analyses:	· · ·			-1

Post-hoc analyses:

<sup>a</sup> Non-Hispanic Blacks are different from Non-Hispanic Whites. <sup>b</sup> Non-Hispanic Whites are different from Non-Hispanic Blacks and Hispanics.

Table 1. Cognivue *Clarity* differences by sociodemographic characteristics in the entire Bio-Hermes study sample (mean [SD] or %).

In the entire study sample, Cognivue Clarity demonstrated discrimination between diagnosis groups (Figure 2) and was able to detect the presence of amyloid (Figure 3).



Figure 2. Mean average score with Cognivue *Clarity* by diagnosis group.



Cognivue *Clarity* was analyzed in 555 biomarker-confirmed participants (297) Figure 4. Mean average score with Cognivue *Clarity* by diagnosis group Controls, 113 MCI, 145 AD) with a mean age of 71.9±6.6 years and 15.6±2.8 years of among participants in the biomarker-confirmed sample. education. 56.4% of participants were female and 21.7% were from underrepresented minorities. Participants had a mean MMSE score of 26.6±3.0, and a mean Cognivue *Clarity* score of 64.0±17.4. Cronbach's alpha was 0.882 (95%CI: 0.867-0.897), Cognivue Clarity was moderately correlated with AD plasma and imaging biomarkers indicating high reliability. (p < .001), while AD biomarkers were moderately correlated with each other (p < .001). Cognivue *Clarity* demonstrated similar associations with AD plasma and imaging As in the entire study sample, Cognivue *Clarity* scores differed between ethnoracial biomarkers as biomarkers have with each other (Table 4).

groups in the biomarker sample (Table 2).

<sup>a</sup> Non-Hispanic Blacks are different from Non-Hispanic Whites.

<sup>b</sup> Hispanics are different from Non-Hispanic Whites. <sup>c</sup> Non-Hispanic Whites are different from Non-Hispanic Blacks and Hispanics.

### Table 2. Cognivue *Clarity* differences by sociodemographic characteristics in the biomarker study sample (mean [SD] or %).

In post-hoc analyses, Cognivue *Clarity* scores were different between the control, MCI-AD, and AD groups, while MMSE scores were different between MCI-AD and AD groups only (Table 3).

	Control	MCI-AD	AD	p-value
Cognivue	71.9 (12.5)	60.1 (16.2)	49.4 (17.4)	<.001
MMSE	28.4 (1.5)	27.0 (1.9)	22.8 (2.6)	<.001
C2N AB42/40	0.101 (.009)	0.089 (.006)	0.091 (.008)	<.001
C2N APS	16.4 (19.4)	59.4 (26.4)	55.9 (29.3)	<.001
Lilly Amyloid PET	0.958 (0.06)	1.376 (0.21)	1.461 (0.21)	<.001
Lilly pTau217	0.171 (0.06)	0.356 (0.23)	0.474 (0.26)	<.001
Quanterix pTau181	16.2 (7.5)	24.1 (11.0)	30.6 (18.5)	<.001
pTau181/AB42	0.327 (0.14)	0.544 (0.27)	0.688 (0.38)	<.001
pTau217/AB42	0.0035 (0.001)	0.0080 (0.005)	0.107 (0.006)	<.001

Table 3. Cognivue *Clarity* differences by biomarker groups (mean [SD]).

Cognivue *Clarity* average and subtest scores robustly differentiated (*p*<0.001) between healthy controls (71.9±12.5), MCI (60.1±16.2) and AD (49.4±17.4) in the

	Cognivue	AB42/40	APS	SUVR	Ptau217	Ptau181	Ptau181/A B42
AB42/40	.253						
APS	339	819					
SUVR	486	455	.584				
Ptau217	426	381	.417	.608			
Ptau181	287	319	.350	.454	.546		
Ptau181/ AB42	307	386	.409	.513	.550	.946	
Ptau217/ AB42	443	411	.442	.629	.958	.503	.588

Table 4. Correlations between Cognivue *Clarity* and biomarkers.

# CONCLUSIONS

Cognivue Clarity discriminated healthy controls from MCI and AD, strongly correlated with AD biomarkers, and had strong psychometric properties. Individuals with Cognivue Clarity average scores <70 were 7.5-fold more likely to be impaired. Cognivue Clarity can be used to screen individuals for cognitive impairment, enrich clinical trial inclusion, and track disease progression.

# FUNDING

This study was funded by grants from AbbVie, Aural Analytics, C2N, Cognivue, EMTherapro, GNS Healthcare, IXICO, Lilly, Linus, Merck, Quanterix, Retispec, Roche, Soma Logic, University of Gothenburg, Global Alzheimer's Platform Foundation<sup>®</sup> (GAP), Alzheimer's Drug Discovery Foundation (ADDF), Gates Ventures, and Pentara Corporation.

# CONTACT

Heather M. Harris. Executive Director, Science, Medical Affairs, and Research, Cognivue Email: hharris@cognivue.com

